
### 5.1

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## Instructions for Use

### DESCRIPTION:

The PleuraFlow® System with FlowGlide™ incorporates a Clearance Apparatus intended to prevent clogging and occlusion of PleuraFlow Chest Tubes used for pleural and mediastinal drainage after cardiothoracic surgery and trauma.

The primary components of the System are the PleuraFlow Chest Tube and the PleuraFlow Clearance Apparatus. (FIG. 1) The PleuraFlow Chest Tube is a silicone chest tube coated with FlowGlide™. It is connected to a Clearance Apparatus, which is connected to the tubing from the drainage canister. (FIG. 3) The PleuraFlow Clearance Apparatus consists of a Guide Tube with a Clearance Wire and Loop that is advanced into the PleuraFlow Chest Tube using a magnetic Shuttle Guide. When indicated, the Clearance Wire and Loop is advanced and retracted within the PleuraFlow Chest Tube to proactively prevent or break up and clear any tube obstructions or clogging to keep the tube patent. Components of the PleuraFlow® System with FlowGlide™ are not made with natural rubber latex.

The device is inserted through the skin adjacent to open surgical incision. The proximal end of the drain is positioned within the operative site prior to repair of the incision. The device’s distal end is attached to an appropriate suction source in order to allow efflux of bloody, serosanguinous, chylous, purulent fluid, and/or other fluids from the operative site that could impair surgical wound healing. The device is indicated for use in cardiothoracic surgical procedures.

#### INDICATIONS:

The PleuraFlow® System with FlowGlide™ is indicated for use during cardiothoracic surgical procedures and chest trauma. Its Active Clearance Technology proactively removes clots formed inside the chest tube to prevent or minimize chest tube occlusion with clot. A patent chest tube enables evacuation of blood and fluid from the operative site after closure of the surgical wound and reduces retained blood.

In the United States the product is indicated for adult and pediatric patients including infant, preadolescent and adolescent patients under clinical settings.

In the European countries the product is indicated for adult and pediatric patients including infants ranging in age from six months old to 18 years old under clinical settings (See Table 1 for additional information).

#### CONTRAINDICATIONS:

The PleuraFlow® System with FlowGlide™ is contraindicated for patients with a history of intolerance to implantable silicone materials.

This device should not be used in proximity to an MRI.

#### SET UP INSTRUCTIONS

##### Prior to Placement:

- Check the package for damage before opening.

#### To open the package:

• Open the Pouch and deliver the contents to the sterile field using aseptic technique.

- Inspect thoroughly, assuring that it is not kinked or otherwise damaged. If there is any damage, replace with a new device.

#### Inserting the PleuraFlow® System with FlowGlide™

• Insert the PleuraFlow Chest Tube into the pleural or mediastinal space according to standard methods.

When a median sternotomy approach is used it is recommended that at least one PleuraFlow® System with FlowGlide™ is used in the anterior mediastinum as the majority of postoperative bleeding occurs in this location.

WARNING: The Chest Tube is not intended for direct contact with the central circulatory system.

Care should be taken to ensure the path of the PleuraFlow Chest Tube is as straight as possible to minimize the resistance of the Clearance Wire and Loop inside the Chest Tube. In some instances, excessive curvature and tortuosity may result in activating the Magnetic Safety Release.

- Secure the PleuraFlow Chest Tube according to standard methods.

• Take care not to constrict the PleuraFlow Chest Tube when securing in place, which may restrict the movement of the Clearance Wire and Loop.

• After insertion, when trimming the PleuraFlow Chest Tube, cut the Chest Tube precisely where indicated by the labeling that indicates "CUT". (FIG. 2)

WARNING: Do not attempt to cut the PleuraFlow Chest Tube shorter than indicated by the "CUT" indicator close to the distal end. Never cut the proximal end (the end having eyelets/drain holes) of the PleuraFlow Chest Tube. This could result in the Clearance Wire and Loop extending beyond the tip of the Chest Tube, which could potentially damage internal structures.

#### Connect the Clearance Apparatus between the PleuraFlow Chest Tube and the drainage tubing (FIG. 3)

• Connect the chest tube to the chest barb, advance tubing all the way onto barb.

• Once the PleuraFlow Clearance Apparatus is connected to the PleuraFlow Chest Tube, advance the external Shuttle Guide toward the proximal barb and Chest Tube. This will advance the Clearance Wire and Loop into the proximal end of the Chest Tube. (FIG. 4)

• Click the Shuttle Guide into the proximal barb housing to park the Clearance Wire and Loop in the proximal end of the Chest Tube.

• Connect the drain barb adapter of the Clearance Apparatus to the drainage tubing that goes to the drainage canister, advance tubing all the way onto barb.

- Secure all connections per hospital protocol.

• Connect the drainage canister to the suction source.

- Maximum vacuum: -40 cmH<sub>2</sub>O

• If a Y-connection is indicated, ensure the Y-junction is placed distal to the Clearance Apparatus. Additional drainage tubing may be used to compensate for length discrepancies.

#### POST INSERTION INSTRUCTIONS

• Confirm tip position of the PleuraFlow Chest Tube according to institution protocol. Although the PleuraFlow Chest Tube material contains a radiopaque stripe to aid in the radiographic visualization of the PleuraFlow Chest Tube, the Clearance Wire and Loop may be left in place to improve radiographic visualization.

#### USE OF THE CLEARANCE APPARATUS

• Only a qualified healthcare practitioner should operate the device.

• When it is indicated to clear the PleuraFlow Chest Tube, the Shuttle Guide is disengaged from the proximal connector, and moved down the Clearance Apparatus, away from the patient and toward the drainage canister tubing. (FIG. 5)

• The Clearance Apparatus should be actuated often in the setting of thick output, such as clotting blood, to ensure the chest tube is patent.

- It is recommend that the device is actuated to clear the PleuraFlow Chest Tube every 15 minutes during the first 8 hours after placement when bleeding is typically more common, then every 30 minutes for the next 16 hours, then every hour thereafter.

• The device should be actuated as needed in addition to these baseline requirements.

• This should be repeated as often as necessary to keep the tube patent and free of any occlusions.

• Each time the Clearance Apparatus is actuated to clear the PleuraFlow Chest Tube; the Clearance Apparatus should be inspected for any clot or occluding material accumulating on the Clearance Wire and Loop.

• If obstructive clot is forming on the Clearance Wire and Loop, steps should be taken to dislodge the clot or fibrinous material stuck to the wire.

• If this cannot be cleared from the wire and is obstructing drainage, the Clearance Wire and Loop should be parked outside the PleuraFlow Chest Tube, in the Clearance Apparatus, by moving the Shuttle Guide to the distal portion of the Clearance Apparatus and leaving it outside of the Chest Tube.

• Traditional methods of chest tube clearance can be carried out at any time, as long as the Clearance Wire and Loop is fully retracted outside of the PleuraFlow Chest Tube.

• When not in use, the Shuttle Guide should be parked by clicking it to the proximal barb, thereby parking the Clearance Wire and Loop in the proximal end of the PleuraFlow Chest Tube. (FIG. 6)

• The Clearance Apparatus should be removed within 5 days or once the bleeding and clotting have ceased, whichever is sooner. This can be done by retracting the Clearance Wire into the Guide Tube and removing the PleuraFlow Chest Tube and the Clearance Apparatus together, if clinically indicated. Alternatively, the Clearance Wire can be retracted into the Guide Tube and the Clearance Apparatus removed, leaving the Chest Tube connected directly to the drainage tubing. The chest tube can then be left in place until removal is clinically indicated up to two weeks from insertion.

#### TROUBLE SHOOTING

- If an obstructive clot appears on the Clearance Wire and Loop, steps should be taken to dislodge the clot into the larger diameter Guide Tube.
- Gently squeeze the wire

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