

PleuraFlow[®] Active Clearance Technology[®] System with **FlowGlide**[™]

The next generation of ClearFlow's Active Clearance Technology

New proprietary treatment of PleuraFlow[®] Active Clearance Technology[®] System with FlowGlide[™] provides a smoother surface inside and outside the chest tube. The improved surface provides:

- **Easier glide path for the clearance apparatus**
- **Tighter bend radius facilitates placement in more tortuous anatomical positions**
- **Smoother chest tube removal may reduce patient discomfort**



Ordering Information

PleuraFlow ACT System with FlowGlide

PleuraFlow ACT System with FlowGlide is available in the same FR sizes and is used in the same manner as the existing PleuraFlow Systems.

CODE	DESCRIPTION	SIZE	NUMBER OF EYELETS	EFFECTIVE DRAINAGE LENGTH (EDL)
PFFG-20	20FR PleuraFlow [®] System with FlowGlide [™]	20 FR	6	4" (10.2 cm)
PFFG-24	24FR PleuraFlow [®] System with FlowGlide [™]	24 FR	6	4" (10.2 cm)
PFFG-28	28FR PleuraFlow [®] System with FlowGlide [™]	28 FR	6	4" (10.2 cm)
PFFG-32	32FR PleuraFlow [®] System with FlowGlide [™]	32 FR	6	4" (10.2 cm)

For Purchase and Re-order of Product

Phone: 1-714-916-5010
Sandra Davalos ext. 814, or Patty Miller ext. 800
Customer Service: 1-714-916-5007

Toll Free: 1-844-257-3569 (1-844-CLR-FLOW)
FAX: 1-714-916-5019
E-Mail orders: customerservice@clearflow.com