


PleuraFlow® XDL

Extended Drainage Length

Enhanced drainage capacity with the confidence of chest tube patency

PleuraFlow® XDL facilitates more effective drainage from harder-to-reach anatomical positions, such as the pleural space and along the diaphragm.



1 **2.5X more drainage length**
than standard devices
(10 inches vs. 4 inches)

2 **2.5X more drainage eyelets**
than standard devices
(15 eyelets vs. 6 eyelets)

3 **Designed to enhance patient comfort**
with our unique reduced friction FlowGlide coating

PLEURAFLOW
ACTIVE CLEARANCE TECHNOLOGY

XDL

The PleuraFlow® ACT® System

is the only FDA cleared device indicated for the removal of retained blood by proactively clearing chest tubes and preventing the retention of blood and fluids in the chest cavity.

CODE	DESCRIPTION	SIZE	NUMBER OF EYELETS	EFFECTIVE DRAINAGE LENGTH (EDL)
PFFG-20 XDL	20FR PleuraFlow® System with FlowGlide® XDL	20FR	15	10" (25.4cm)
PFFG-24 XDL	24FR PleuraFlow® System with FlowGlide® XDL	24FR	15	10" (25.4cm)

 **CLEARFLOW**
CLEARING THE PATH TO RECOVERY

Copyright © 2018 ClearFlow, Inc.
All Rights Reserved.

™ Trademarks and * Registered Trademarks of ClearFlow, Inc.

Manufactured For:
ClearFlow, Inc.
1630 S. Sunkist St., Suite E
Anaheim, CA, US, 92806
www.clearflow.com

Contact Customer Service:
714-916-5007
customerservice@clearflow.com

ML239-A